



**RFID SECURITY**  
**INCREASE SAFETY**  
**FOR YOUR PATIENTS**  
**AND RESIDENTS**

A photograph of an elderly woman with dark hair pulled back, wearing a light green cardigan over a white collared shirt. She is standing and holding the handles of a silver metal walker. The background is a bright, indoor setting, possibly a hospital or care facility, with a window and blinds visible. The image is framed by teal geometric shapes in the corners.

# **Patient Wander Alert System**

Protrac iD Patient Wandering Alert Systems detect wandering patients wearing Protrac iD RFID tags when they approach monitored doors or areas both inside and outside.

Depending on the settings and the mode of operation selected the PWAS triggers alarms, locks or open doors, alerts staff and communicates with other systems such as Nurse Call, Access Control, Radio Paging, DECT Phones, Protrac iD AssetTrac or BabyGuard™ Software.

Protrac iD Patient Wandering Alert System is available as a 2R5000 or 2R5500.

Protrac iD tags are individually identified and transmit a "heart beat" radio signal continuously with a battery life exceeding 3 years.

### ***Benefits of the Protrac iD Patient Wander Alert System:***


- Give more freedom to patients/residents
- Integrated into most new & existing Nurse Call Systems
- Indoor & Outdoor alert monitoring
- Up to 150m detection range
- No ugly portals
- Wireless
- Door lock
- Staff override



**2R5000PW**



**2R5500PW**



**1RN350 Patch Antenna**



**2T7000 Bracelet Tag**

Distributed by:



Designed and Manufactured in Australia.

All rights reserved. Specifications subject to change without notice.

[www.protracid.com.au](http://www.protracid.com.au)